

IHE Integration Statement			
Vendor	Product Name	Version	Date
FUJIFILM Corporation	Flash IIP (CR-IR346/348CL)	Flex UI v6.0	4/16/2010
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF)	Acquisition Modality	Broad Worklist Query Modality Group Case	
Patient Information Reconciliation (PIR)			
Audit Trail & Node Authentication (ATNA)	Secure Node	N/A	
Portable Data for Imaging (PDI)	Portable Media Creator	-	
Internet address for vendor's IHE information: http://www.fujimed.com/sub/product_bulletins.html			
Links to Standards Conformance Statements for the Implementation			
HL7			
DICOM	http://www.fujimed.com/sub/product_bulletins.html		
Links to general information on IHE			
In North America: www.ihe.net	In Europe: www.ihe-europe.net	In Japan: www.jira-net.or.jp/ihe-j	